


## Dentrix Trojan CD Update Traditional Database Update Instructions


1. Open the **Dentrix Trojan Update Utility** If you do not have an icon for this on your desktop, Click on **Start, Programs, Dentrix Trojan Viewer 6.0, Dentrix Trojan Update Utility**



You may have an icon like this already on your desktop.

2. Click on **Database Update**
3. *You will see a screen alerting you that the database will be modified, Simply click <OK>.*

4. Select the Insurance Plans and choose the location of the updates. In the **Install Update From:** box, Type **D:>** for your CD drive (your drive letter may vary), and click <OK>.



5. Dentrix will now start to process the update.
6. Once the update is complete you will see a box that will say **DONE** Click <OK>